

EBP06173 - IMMUNIZING PEPTIDE

SKU: EBP06173



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_008970

Immunogen Peptide with sequence C-RVLIQFLKSSRLI, from the C Terminus of the protein sequence according to NP_008970.

Product Comments The peptide was used in the production of EB06173 - Goat Anti-PTPD1 / PTPN21 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-RVLIQFLKSSRLI

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png