

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06182 - IMMUNIZING PEPTIDE

SKU: EBP06182

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP_003244; NP_001340613.1; NP_001340614.1; NP_001340615.1; NP_001340616.1 ID

Peptide with sequence C-EDFAPSRKLNTEI, from the C Terminus of the protein sequence according to Immunogen

NP 003244; NP 001340613.1; NP 001340614.1; NP 001340615.1; NP 001340616.1.

The peptide was used in the production of EB06182 - Goat Anti-TIAM1 Antibody. It is currently untested in **Product**

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Comments blocking applications.

Peptide

C-EDFAPSRKLNTEI Sequence

Shipping

Product

Instructions Ambient

Grade