

EBP06196 - IMMUNIZING PEPTIDE

SKU: EBP06196



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005707.1; NP_974197.1; NP_974199.1; NP_974196.1; NP_974198.1; NP_974223.1
Immunogen	Peptide with sequence C-GAIGDAKVGRY, from the C Terminus of the protein sequence according to NP_005707.1; NP_974197.1; NP_974199.1; NP_974196.1; NP_974198.1; NP_974223.1.
Product Comments	The peptide was used in the production of EB06196 - Goat Anti-GIPC1 / NIP Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GAIGDAKVGRY
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)