

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06201 - IMMUNIZING PEPTIDE**

**SKU:** EBP06201

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** NP\_079120.1; NP\_001290205.1; NP\_001317401.1 ID

Immunogen Peptide with sequence C-ARAIIREGSLEGS, from the C Terminus of the protein sequence according to NP\_079120.1; NP\_001290205.1; NP\_001317401.1.

The peptide was used in the production of EB06201 - Goat Anti-Pericentrin 1 / NUP85 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-ARAIIREGSLEGS Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade