

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06202 - IMMUNIZING PEPTIDE

SKU: EBP06202

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_006836.1 ID

 $\begin{tabular}{ll} \textbf{Immunogen} \\ \textbf{NP_006836.1}. \\ \end{tabular} Peptide with sequence AMDSSLQARLFPG-C, from the N Terminus of the protein sequence according to NP_006836.1. \\ \end{tabular}$

The peptide was used in the production of EB06202 - Goat Anti-MCAK / KIF2C Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

AMDSSLQARLFPG-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade