

EBP06207 - IMMUNIZING PEPTIDE

SKU: EBP06207



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_002690.3

Immunogen Peptide with sequence ATTAQYLPRGPGG-C, from the N Terminus of the protein sequence according to NP_002690.3.

Product Comments The peptide was used in the production of EB06207 - Goat Anti-OCT6 / SCIP / POU3F1 Antibody. It is currently untested in blocking applications.

Peptide Sequence ATTAQYLPRGPGG-C

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png