

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06213 - IMMUNIZING PEPTIDE

SKU: EBP06213

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 036478.2; NP 057637.2; NP 714940.1; NP 714941.1 ID

Peptide with sequence C-RKEQERSFRITFD, from the C Terminus of the protein sequence according to Immunogen

NP 036478.2; NP 057637.2; NP 714940.1; NP 714941.1.

The peptide was used in the production of EB06213 - Goat Anti-NUP62 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide C-RKEQERSFRITFD Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade