

EBP06213 - IMMUNIZING PEPTIDE

SKU: EBP06213



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_036478.2; NP_057637.2; NP_714940.1; NP_714941.1
Immunogen	Peptide with sequence C-RKEQERSFRITFD, from the C Terminus of the protein sequence according to NP_036478.2; NP_057637.2; NP_714940.1; NP_714941.1.
Product	The peptide was used in the production of EB06213 - Goat Anti-NUP62 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-RKEQERSFRITFD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)