

EBP06217-B - IMMUNIZING PEPTIDE

SKU: EBP06217-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_055241,1; NP_001269383,1
Immunogen	Peptide with sequence SLLQSCIDLFKNN,, from the C Terminus of the protein sequence according to NP_055241,1; NP_001269383,1,
Product	The peptide was used in the production of EB06217-B - Goat Anti-SNX5, Biotinylated Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	SLLQSCIDLFKNN,
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)