

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06217-B - IMMUNIZING PEPTIDE**

**SKU:** EBP06217-B

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 055241,1; NP 001269383,1 ID

Peptide with sequence SLLQSCIDLFKNN,, from the C Terminus of the protein sequence according to Immunogen

NP\_055241,1; NP\_001269383,1,

The peptide was used in the production of EB06217-B - Goat Anti-SNX5, Biotinylated Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

SLLQSCIDLFKNN, Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade