

EBP06220 - IMMUNIZING PEPTIDE

SKU: EBP06220



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_057308

Immunogen Peptide with sequence C-EKLRQALSRFPVM, from the C Terminus of the protein sequence according to NP_057308.

Product Comments The peptide was used in the production of EBP06220 - Goat Anti-SNX9 / SH3PX1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-EKLRQALSRFPVM

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png