

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06223 - IMMUNIZING PEPTIDE

SKU: EBP06223

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

ID

Accession NP_071919.2; NP_001243479.1

 $\label{eq:local_$

Product The peptide was used in the production of EB06223 - Goat Anti-PALS / MPP5 Antibody. It is currently untested in

Comments blocking applications.

Peptide

NKASPPFPLISNAQ-C Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade