

EBP06233 - IMMUNIZING PEPTIDE

SKU: EBP06233

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.**Unit Size** 100 µg**Storage Instructions** Shipped at ambient temperature, store at -20°C.**Accession ID** NP_002257.1**Immunogen** Peptide with sequence C-QVQDGAPGTFNF, from the C Terminus of the protein sequence according to NP_002257.1.**Product Comments** The peptide was used in the production of EB06233 - Goat Anti-KPNA2 / IPOA1 Antibody. It is currently untested in blocking applications.**Peptide Sequence** C-QVQDGAPGTFNF**Shipping Instructions** Ambient**Product Grade** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)