

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06242 - IMMUNIZING PEPTIDE

SKU: EBP06242

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_002262 ID

Immunogen Peptide with sequence PEDQVGKLEATEN-C, from the N Terminus of the protein sequence according to NP_002262. Product The peptide was used in the production of EB06242 - Goat Anti-KPNB3 / RANBP5 Antibody. It is currently untested

Comments in blocking applications.

Peptide

Sequence

PEDQVGKLEATEN-C

Shipping Ambient Instructions

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png