

# EBP06245 - IMMUNIZING PEPTIDE

**SKU:** EBP06245



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_005692

**Immunogen** Peptide with sequence SHSPDHPGCLMEN, from the C Terminus of the protein sequence according to NP\_005692.

**Product Comments** The peptide was used in the production of EBP06245 - Goat Anti-Snurportin1 Antibody. It is currently untested in blocking applications.

**Peptide Sequence** SHSPDHPGCLMEN

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)