

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06248 - IMMUNIZING PEPTIDE

SKU: EBP06248

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_149094.3 ID

The peptide was used in the production of EB06248 - Goat Anti-Rhophilin 2 (N Terminus) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

TDALLPAAPQPLEK-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade