

EBP06258 - IMMUNIZING PEPTIDE

SKU: EBP06258



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_075067
Immunogen	Peptide with sequence C-NGTSWIEALMEKPF, from the C Terminus of the protein sequence according to NP_075067.
Product Comments	The peptide was used in the production of EB06258 - Goat Anti-VPS33A Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NGTSWIEALMEKPF
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png