



# EBP06260 - IMMUNIZING PEPTIDE

**SKU:** EBP06260



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage</b>	
<b>Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_694985.28; NP_006738.3
<b>Immunogen</b>	Peptide with sequence C-ASKQLGSPTADLA, from the C Terminus of the protein sequence according to NP_694985.28; NP_006738.3.
<b>Product</b>	The peptide was used in the production of EB06260 - Goat Anti-SPA1 / SIPA1 Antibody. It is currently untested
<b>Comments</b>	in blocking applications.
<b>Peptide Sequence</b>	C-ASKQLGSPTADLA
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)