

# EBP06266 - IMMUNIZING PEPTIDE

**SKU:** EBP06266



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_000174.1

**Immunogen** Peptide with sequence TILTYPFKNLPT-C, from the N Terminus of the protein sequence according to NP\_000174.1.

**Product Comments** The peptide was used in the production of EB06266 - Goat Anti-HADHB Antibody. It is currently untested in blocking applications.

**Peptide Sequence** TILTYPFKNLPT-C

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)

## DOCUMENTS

- [Data Sheet](#)