

EBP06271 - IMMUNIZING PEPTIDE

SKU: EBP06271



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_061830.3

Immunogen Peptide with sequence C-PQPRPRSLASETN, from the C Terminus of the protein sequence according to NP_061830.3.

Product Comments The peptide was used in the production of EB06271 - Goat Anti-SH3BP1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-PQPRPRSLASETN

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)