

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06277 - IMMUNIZING PEPTIDE

SKU: EBP06277

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_110379.2; NP_001008897.1

 $\begin{tabular}{ll} \textbf{Immunogen} \\ \textbf{NP_110379.2}; \ \textbf{NP_001008897.1}. \\ \end{tabular}$

The peptide was used in the production of EB06277 - Goat Anti-TCP1 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-SYEDAVHSGALND Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade