

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06284 - IMMUNIZING PEPTIDE

SKU: EBP06284

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001035948.1; NP_036517 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \frac{\text{Peptide with sequence ATFPPATSAPQQP-C, from the N Terminus of the protein sequence according to NP_001035948.1; NP_036517.}$

The peptide was used in the production of EB06284 - Goat Anti-P8 / NUPR1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

ATFPPATSAPQQP-C Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade