

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06304 - IMMUNIZING PEPTIDE

SKU: EBP06304

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_003783.1; \ NP_694880.1; \ NP_001268226.1; \ NP_001268227.1; \ NP_001268228.1$ ID

Peptide with sequence AESDWDTVTVLRK-C, from the N Terminus of the protein sequence according to

Immunogen NP_003783.1; NP_694880.1; NP_001268226.1; NP_001268227.1; NP_001268228.1.

The peptide was used in the production of EB06304 - Goat Anti-EDF1 / MBF1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

AESDWDTVTVLRK-C Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade