



EBP06304 - IMMUNIZING PEPTIDE

SKU: EBP06304



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003783.1; NP_694880.1; NP_001268226.1; NP_001268227.1; NP_001268228.1
Immunogen	Peptide with sequence AESDWDTVTVLRK-C, from the N Terminus of the protein sequence according to NP_003783.1; NP_694880.1; NP_001268226.1; NP_001268227.1; NP_001268228.1.
Product Comments	The peptide was used in the production of EB06304 - Goat Anti-EDF1 / MBF1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	AESDWDTVTVLRK-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png