

EBP06321 - IMMUNIZING PEPTIDE

SKU: EBP06321



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_055126.1; NP_001230718.1; NP_001338187.1; NP_001338190.1
Immunogen	Peptide with sequence PQAFLLSIHEPA-C, from the N Terminus of the protein sequence according to NP_055126.1; NP_001230718.1; NP_001338187.1; NP_001338190.1.
Product	The peptide was used in the production of EB06321 - Goat Anti-SMUG1 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	PQAFLLSIHEPA-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)