

EBP06329 - IMMUNIZING PEPTIDE

SKU: EBP06329



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_036354.1; NP_001041636.1; NP_001041637.1; NP_001041638.1; NP_001121897.1
Immunogen	Peptide with sequence C-HISTDAHSLNSAAQ, from the C Terminus of the protein sequence according to NP_036354.1; NP_001041636.1; NP_001041637.1; NP_001041638.1; NP_001121897.1.
Product Comments	The peptide was used in the production of EBP06329 - Goat Anti-MUTYH Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-HISTDAHSLNSAAQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)