

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06329 - IMMUNIZING PEPTIDE

SKU: EBP06329

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 036354.1; NP 001041636.1; NP 001041637.1; NP 001041638.1; NP 001121897.1 ID

Peptide with sequence C-HISTDAHSLNSAAQ, from the C Terminus of the protein sequence according to Immunogen

NP_036354.1; NP_001041636.1; NP_001041637.1; NP_001041638.1; NP_001121897.1.

The peptide was used in the production of EB06329 - Goat Anti-MUTYH Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide C-HISTDAHSLNSAAQ Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade