

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06331 - IMMUNIZING PEPTIDE

SKU: EBP06331

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_005093.1 ID

 $\begin{tabular}{ll} \textbf{Immunogen} & Peptide with sequence GERRGRDGLGAPT-C, from the N Terminus of the protein sequence according to NP_005093.1. \end{tabular}$

The peptide was used in the production of EB06331 - Goat Anti-FEZ2 / Zygin 2 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

GERRGRDGLGAPT-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade