

EBP06334 - IMMUNIZING PEPTIDE

SKU: EBP06334



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_112497.2; NP_001159438.1
Immunogen	Peptide with sequence C-SWARGPPKFESVP, from the C Terminus of the protein sequence according to NP_112497.2; NP_001159438.1.
Product	The peptide was used in the production of EB06334 - Goat Anti-NIR1 / PITPNM3 Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-SWARGPPKFESVP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)