

EBP06339-B - IMMUNIZING PEPTIDE

SKU: EBP06339-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001092692.1; NP_001092693.1; NP_079379.2; AAK77221.1
Immunogen	Peptide with sequence C-KTIVDALHSTSGN., from the C Terminus of the protein sequence according to NP_001092692.1; NP_001092693.1; NP_079379.2; AAK77221.1.
Product	The peptide was used in the production of EBP06339-B - Goat Anti-CMG1 / CCDC2 / IFT74, Biotinylated Antibody.
Comments	It is currently untested in blocking applications.
Peptide Sequence	C-KTIVDALHSTSGN.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)