

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06339-B - IMMUNIZING PEPTIDE

SKU: EBP06339-B

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP_001092692.1; NP_001092693.1; NP_079379.2; AAK77221.1 ID

Peptide with sequence C-KTIVDALHSTSGN., from the C Terminus of the protein sequence according to Immunogen

NP 001092692.1; NP 001092693.1; NP 079379.2; AAK77221.1.

The peptide was used in the production of EB06339-B - Goat Anti-CMG1 / CCDC2 / IFT74, Biotinylated Antibody. **Product**

Comments It is currently untested in blocking applications.

Peptide

C-KTIVDALHSTSGN. Sequence

Shipping

Instructions Ambient

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png