

# EBP06341 - IMMUNIZING PEPTIDE

**SKU:** EBP06341



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_065881.1

**Immunogen** Peptide with sequence VVGCSLATGFSPH, from the C Terminus of the protein sequence according to NP\_065881.1.

**Product Comments** The peptide was used in the production of EBP06341 - Goat Anti-FENS1 / WDFY1 Antibody. It is currently untested in blocking applications.

**Peptide Sequence** VVGCSLATGFSPH

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)

## DOCUMENTS

- [Data Sheet](#)