

EBP06342 - IMMUNIZING PEPTIDE

SKU: EBP06342



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_002430

Immunogen Peptide with sequence C-EFNMEETQTSLH, from the C Terminus of the protein sequence according to NP_002430.

Product Comments The peptide was used in the production of EBP06342 - Goat Anti-MSH3 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-EFNMEETQTSLH

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png