

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06343 - IMMUNIZING PEPTIDE

SKU: EBP06343

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

Accession NP_852122.1; NP_000429.2; NP_039230.1; NP_852123.1; NP_852124.1; NP_852126.1; NP_852125.1; ID NP_001120838.1

Peptide with sequence TTLAGAVPRMM-C, from the N Terminus of the protein sequence according to

Immunogen NP 852122.1; NP 000429.2; NP 039230.1; NP 852123.1; NP 852124.1; NP 852126.1; NP 852125.1;

NP 001120838.1.

The peptide was used in the production of EB06343 - Goat Anti-PAX3 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

TTLAGAVPRMM-C Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade