

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06344 - IMMUNIZING PEPTIDE

SKU: EBP06344

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_004292.1; NP_817126.1; NP_057272.1 ID

Immunogen Peptide with sequence YEEGGKQCVERKCP, from the C Terminus of the protein sequence according to NP_004292.1; NP_817126.1; NP_057272.1.

The peptide was used in the production of EB06344 - Goat Anti-BAF53A and BAF53B Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

YEEGGKQCVERKCP Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade