

EBP06344 - IMMUNIZING PEPTIDE

SKU: EBP06344



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004292.1; NP_817126.1; NP_057272.1
Immunogen	Peptide with sequence YEEGGKQCVERKCP, from the C Terminus of the protein sequence according to NP_004292.1; NP_817126.1; NP_057272.1.
Product Comments	The peptide was used in the production of EB06344 - Goat Anti-BAF53A and BAF53B Antibody. It is currently untested in blocking applications.
Peptide Sequence	YEEGGKQCVERKCP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)