

EBP06351 - IMMUNIZING PEPTIDE

SKU: EBP06351



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002432.1; NP_079535.4; NP_751897.1
Immunogen	Peptide with sequence C-MSQEVLPAAATIL, from the C Terminus of the protein sequence according to NP_002432.1; NP_079535.4; NP_751897.1.
Product Comments	The peptide was used in the production of EBP06351 - Goat Anti-MSH5 (C-term.) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-MSQEVLPAAATIL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)