

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06351 - IMMUNIZING PEPTIDE

SKU: EBP06351

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_002432.1; NP_079535.4; NP_751897.1 ID

The peptide was used in the production of EB06351 - Goat Anti-MSH5 (C-term.) Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

Sequence

C-MSQEVLPAATSIL

Shipping

Ambient Instructions

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png