

EBP06352 - IMMUNIZING PEPTIDE

SKU: EBP06352



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_003060; NP_620604

Immunogen Peptide with sequence C-ATKTPMVKFSAFS, from the C Terminus of the protein sequence according to NP_003060; NP_620604.

Product Comments The peptide was used in the production of EB06352 - Goat Anti-SNF2L1 / SMARCA1 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-ATKTPMVKFSAFS

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png