

EBP06353 - IMMUNIZING PEPTIDE

SKU: EBP06353



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001122321.1; NP_003063.2; NP_001122317.1; NP_001122318.1; NP_001122319.1; NP_001122320.1
Immunogen	Peptide with sequence STDPPLGGTPRP-C, from the N Terminus of the protein sequence according to NP_001122321.1; NP_003063.2; NP_001122317.1; NP_001122318.1; NP_001122319.1; NP_001122320.1.
Product Comments	The peptide was used in the production of EB06353 - Goat Anti-BRG1 / SMARCA4 Antibody. It is currently untested in blocking applications.
Peptide Sequence	STDPPLGGTPRP-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)