

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06353 - IMMUNIZING PEPTIDE

SKU: EBP06353

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

ID

 $NP_001122321.1; \ NP_003063.2; \ NP_001122317.1; \ NP_001122318.1; \ NP_001122319.1; \ NP_001122320.1 \ NP_001122319.1; \ NP_001122319.1$

Peptide with sequence STPDPPLGGTPRP-C, from the N Terminus of the protein sequence according to

Immunogen NP_001122321.1; NP_003063.2; NP_001122317.1; NP_001122318.1; NP_001122319.1; NP_001122320.1.

The peptide was used in the production of EB06353 - Goat Anti-BRG1 / SMARCA4 Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

Sequence

STPDPPLGGTPRP-C

Shipping

Ambient Instructions

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$