



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06376 - IMMUNIZING PEPTIDE

SKU: EBP06376

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_473380.1 ID

Immunogen Peptide with sequence KRSQRPNKCSNP, from the C Terminus of the protein sequence according to NP_473380.1.

Product The peptide was used in the production of EB06376 - Goat Anti-FOXP3 / Scurfin (Mouse) Antibody. It is currently

Comments untested in blocking applications.

Peptide

KRSQRPNKCSNP Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade