

EBP06377 - IMMUNIZING PEPTIDE

SKU: EBP06377



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_002037.2

Immunogen Peptide with sequence C-HQVVSSDFNSDT, from the C Terminus of the protein sequence according to NP_002037.2.

Product Comments The peptide was used in the production of EB06377 - Goat Anti-GAPDH (C Terminus) Loading Control Antibody. It is currently untested in blocking applications.

Peptide Sequence C-HQVVSSDFNSDT

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)