



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06394 - IMMUNIZING PEPTIDE

SKU: EBP06394

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_001620.2; NP_001191335.1 ID

Immunogen Peptide with sequence C-KTISTTISPLLLIP, from the C Terminus of the protein sequence according to NP_001620.2; NP_001191335.1.

The peptide was used in the production of EB06394 - Goat Anti-FLAP / ALOX5AP (C Term) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-KTISTTISPLLLIP Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade