

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06435 - IMMUNIZING PEPTIDE**

**SKU:** EBP06435

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_065144; CAC38108 ID

Immunogen Peptide with sequence CPEAGEEQPQV, from the C Terminus of the protein sequence according to NP\_065144; CAC38108.

The peptide was used in the production of EB06435 - Goat Anti-XAGE1 (Isoforms 1a, 1b and 1c) Antibody. It is Product

**Comments** currently untested in blocking applications.

**Peptide** 

**CPEAGEEQPQV** Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade