

EBP06436 - IMMUNIZING PEPTIDE

SKU: EBP06436



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_872584.2; NP_079043.3
Immunogen	Peptide with sequence SVAIRKRSWEEHC, from the N Terminus of the protein sequence according to NP_872584.2; NP_079043.3.
Product Comments	The peptide was used in the production of EB06436 - Goat Anti-DLC1 (Isoforms 1 and 3) Antibody. It is currently untested in blocking applications.
Peptide Sequence	SVAIRKRSWEEHC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)