

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06442 - IMMUNIZING PEPTIDE

SKU: EBP06442

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_006648.1; NP_001307341.1 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \begin{array}{l} \text{Peptide with sequence SQLVECVPNFSEGKNQ, from the N Terminus of the protein sequence according to NP_006648.1; NP_001307341.1.} \end{array}$

The peptide was used in the production of EB06442 - Goat Anti-58K Golgi protein(N-Term)/FTCD Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

SQLVECVPNFSEGKNQ Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade