

EBP06443 - IMMUNIZING PEPTIDE

SKU: EBP06443



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006648.1; NP_996848.1; NP_001307341.1
Immunogen	Peptide with sequence CLREQGRGKDQPGRL, from the internal region of the protein sequence according to NP_006648.1; NP_996848.1; NP_001307341.1.
Product Comments	The peptide was used in the production of EB06443 - Goat Anti-58KGolgi protein(Internal)/FTCD Antibody. It is currently untested in blocking applications.
Peptide Sequence	CLREQGRGKDQPGRL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)