

EBP06468 - IMMUNIZING PEPTIDE

SKU: EBP06468



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002934.1; NP_599022.1; NP_599023.1; NP_599024.1
Immunogen	Peptide with sequence CTSEFEPAMQIDG, from the C Terminus of the protein sequence according to NP_002934.1; NP_599022.1; NP_599023.1; NP_599024.1.
Product	The peptide was used in the production of EBP06468 - Goat Anti-RORA Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	CTSEFEPAMQIDG
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png