



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06476 - IMMUNIZING PEPTIDE**

**SKU:** EBP06476

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_853523.1 ID

Immunogen Peptide with sequence YFPTPEGHSSTEAC, from the C Terminus of the protein sequence according to NP\_853523.1.

Product The peptide was used in the production of EB06476 - Goat Anti-Schlafen 8 (Mouse) Antibody. It is currently

**Comments** untested in blocking applications.

**Peptide** 

Sequence

YFPTPEGHSSTEAC

**Shipping** 

Ambient Instructions

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png